

Frequently Asked Questions

PATIENT TYPE

Can the Sensium patch be used on maternity patients?

Yes, the patch can be used on maternity patients following the instructions in the Patch IFU

Can the Sensium patch be used on children?

The Patch has not yet received regulatory approval for use on children under 18.

COMMUNICATION

Can the Sensium patch communicate through WiFi?

The system can communicate via Wi-Fi provided that a robust and reliable Wi-Fi connection is available.

Our WiFi and IT is poor, can we still use it?

Please refer to our minimum technical requirements document to see if your hospital IT infrastructure is sufficient

Can the Sensium system communicate through Bluetooth?

No, the patch does not operate using Bluetooth.

THE PATCH

How long does the Sensium patch last?

The Patch has a battery life of 5 days

Can we change the battery and reuse the patch?

The patch is designed for single-patient use to avoid the risk of cross-infection

Can a patient shower whilst wearing the patch?

The Patch is showerproof; however, we recommend that you change the electrodes after the shower.

COSTS & INSTALLATION

How long does the Sensium system take to install?

Installation time will vary depending on the size of the hospital/ward and support from internal IT, but can be as short as 4 weeks

How long does the Sensium system cost?

The Patch, Bridge, licensing and installation costs will vary depending on the size of the installation and type of deployment